



8/32" MOUNTING HOLES

MOTOR DRIVER CABLE  
PART NO. I-22640-CE  
3 FT LONG

POWER CABLE 24VDC

MOUNTING BRACKET  
PART NO. I-22766

MANUAL CONTROLLER  
PART NO. I-22679

**DIMENSIONS:**  
58 x 92 x 27 mm  
2.27 x 3.62 x 1.06 inch

REMOTE CONTROL CABLE  
PART NO. I-22696  
6 FT LONG

MANUAL CONTROLLER CABLE  
PART NO. I-22696-01  
3 FT LONG

**VICI** Valco Instruments Co. Inc.

MULTIPOSITION ACTUATOR  
CONTROL MODULE

CONTROL MODULE  
BOTTOM SIDE

CONTROL MODULE  
PART NO. FOR EMH (High speed): EMHCA-CE  
PART NO. FOR EMT (High torque): EMTCA-CE

**DIMENSIONS:**  
94 x 114 x 38.1 mm  
3.7 x 4.5 x 1.5 inch

Dimension in mm [Inch]

<small>THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.</small>		Model: <b>MTEFAH</b>	Dimensional tolerances: GENERAL TOLERANCES
Änderung: Modifikation:	a) 25.03.13 FP #716, chg Mon. controller design	Projektion: PROJEKTION	Alternative: Toleranzwerte
PART NO.:		Zeichnung: DRAWN	Date: 10.09.1998 F.P.
		Maßstab: SCALE: 1:5.1	Zeichnung: DRAWN: 25.03.2013 FCo
		Material: MAT:	Zeichnung: DRAWN: 506019.DWG
		Material: Star code:	Zeichnung: DRAWN: KUNDE
		Zeichnung: DRAWN:	Zeichnung: DRAWN:
		<b>VICI</b>	Vici AG International, CH-6214 Scherikon <b>1-GE-5060.19 a</b>