



PN#: Material:  
\*PK PEEK




SHOWN ZN1PK

Part#	Description	L		X	
		mm	Inch	mm	Inch
ZN1PK	Nut, ZDV: 1/16", standard	11.43	0.450"	1/4" Hex	
MZN1PK	Nut, ZDV: 1/16", long	15.88	0.625"	1/4" Hex	
LZN1PK	Nut, ZDV: 1/16", long	22.23	0.875"	1/4" Hex	
ZN2PK	Nut, ZDV: 1/8", standard	15.88	0.625"	3/8" Hex	



PN#: Material:  
\*FPK PEEK  
\*FKF CTFE  
\*FD DELRIN



SHOWN ZN1FPK

Part#	Description	L		X	
		mm	Inch	mm	Inch
ZN.5FPK	Nut, ZDV: 1/32", standard	10.67	0.420"	∅6.35	∅0.250"
LZN.5FPK	Nut, ZDV: 1/32", long	13.72	0.540"	∅6.35	∅0.250"
ZN1FPK	Nut, ZDV: 1/16", standard	22.23	0.875"	∅9.65	∅0.380"
ZN1FKF	Nut, ZDV: 1/16", standard	22.23	0.875"	∅9.65	∅0.380"
ZN1FD	Nut, ZDV: 1/16", standard	22.23	0.875"	∅9.65	∅0.380"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material  
MATERIAL


Dimensionen  
DIMENSIONS

Aenderungen:  
MODIFICATIONS:

Projektion  
PROJECTION

Allgemein Toleranzen  
GENERAL TOLERANCES

PART NO.:



SN EN 22768-1 / ISO 2768

# Polymeric nuts High pressure

General dimensions

Masstab  
SCALE

Gezeichnet  
DRAWN

15.12.2005

F.P.

2:1

Geprüft / Freigabe  
CHECKED/ RELEASE

15.12.2005

FCo

CAD Nr.  
CAD NO.

506058.idw

Unterverz.:  
SUB DIREC.:

KUNDE

Zeichnungs Nr. / DRAWING NO.

4 - GE - 5060.58



VICI AG International, CH-6214 Schenkon