



Shown CFPA2PK

Shown CEPA2PK

Shown CPA2PK

Dimension in mm [Inch]

Part#	Description	Material	L	L1	X	Thread
CPA2PK	Female 1/4-28 to male NPT	PEEK	20.32 [.80"]	11.4 [.45"]	Ø15.9 [Ø5/8"]	1/8-27 NPT
CPA2KF	Female 1/4-28 to male NPT	CTFE	20.32 [.80"]	11.4 [.45"]	Ø15.9 [Ø5/8"]	1/8-27 NPT
CPA4PK	Female 1/4-28 to male NPT	PEEK	23.11 [.91"]	14.2 [.56"]	Ø19.05 [Ø3/4"]	1/4-18 NPT
CPA4KF	Female 1/4-28 to male NPT	CTFE	23.11 [.91"]	14.2 [.56"]	Ø19.05 [Ø3/4"]	1/4-18 NPT
CEPA2PK	Male 1/4-28 to male NPT	PEEK	28.7 [1.13"]	11.4 [.45"]	11.10 [.437"] SQ.	1/8-27 NPT
CEPA2KF	Male 1/4-28 to male NPT	CTFE	28.7 [1.13"]	11.4 [.45"]	11.10 [.437"] SQ.	1/8-27 NPT
CEPA4PK	Male 1/4-28 to male NPT	PEEK	31.0 [1.22"]	14.2 [.56"]	14.27 [.562"] SQ.	1/4-18 NPT
CEPA4KF	Male 1/4-28 to male NPT	CTFE	31.0 [1.22"]	14.2 [.56"]	14.27 [.562"] SQ.	1/4-18 NPT
CFPA2PK	Female 1/4-28 to female NPT	PEEK	19.56 [.77"]	9.14 [.36"]	Ø15.9 [Ø5/8"]	1/8-27 NPT
CFPA2KF	Female 1/4-28 to female NPT	CTFE	19.56 [.77"]	9.14 [.36"]	Ø15.9 [Ø5/8"]	1/8-27 NPT
CFPA4PK	Female 1/4-28 to female NPT	PEEK	24.89 [.98"]	14.5 [.57"]	Ø19.05 [Ø3/4"]	1/4-18 NPT
CFPA4KF	Female 1/4-28 to female NPT	CTFE	24.89 [.98"]	14.5 [.57"]	Ø19.05 [Ø3/4"]	1/4-18 NPT

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: a) 10.05.16 FP #998, changed NPT thread length ECN# 21261  
 MODIFICATIONS:

PART NO.:

**Pipe adapters  
1/4-28 UNF to NPT**

General dimensions

**VICI** VICI AG International, CH-6214 Schenkon

Masstab SCALE	Gezeichnet DRAWN	04.06.2009	F.P.
2:1	Geprüft / Freigabe CHECKED/ RELEASE	10.05.2016	FCo
	CAD Nr. CAD NO.	506350.idw	
Unterverz.: SUB DIREC.:	KUNDE		
Zeichnungs Nr. / DRAWING NO.	3 - GE - 5063.50 a		

Material MATERIAL      Dimensionen DIMENSIONS  
 Projektion PROJECTION      Allgemein Toleranzen GENERAL TOLERANCES  
 SN EN 22768-1/ ISO 2768