



Dimension in mm [Inch]

Prod No	Description	Bore	L	L1
ASFVHT	Air actuated with 2" standoff	0.50 mm	~ 130.5 [~ 5.14 "]	50.55 [1.99"]
ASFVLHT	Air actuated with 2" standoff	0.75 mm	~ 130.5 [~ 5.14 "]	50.55 [1.99"]
ASFVHT4	Air actuated with 4" standoff	0.50 mm	~ 181.5 [~ 7.14 "]	101.35 [3.99"]
ASFVLHT4	Air actuated with 4" standoff	0.75 mm	~ 181.5 [~ 7.14 "]	101.35 [3.99"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: a) 22.05.12 FP #689, changed description to 2" standoff
MODIFICATIONS:

PART NO.: ASFV*

Gezeichnet: 24.02.2010 F.P.
Geprüft / Freigabe: 22.05.2012 FCo
CHECKED/RELEASE
CAD Nr.: 506359.idw
CAD NO.

Material: DIMENSIONEN
Allgemein Toleranzen: GENERAL TOLERANCES
SN EN 22768-1/ ISO 2768

Projection: PROJEKTION

Unterverz.: KUNDE
SUB DIREC.: Zeichnungs Nr. / DRAWING NO.

Prime / purge valve, 1/16", air actuated, high temperature

General dimensions

Masstab: 1:1
SCALE

3 - GE - 5063.59 a

VICI® VICI AG International, CH-6214 Schenkon