



Part#	Sample size	Barrel OD. X1	Barrel ID. X2	Front/ rear sleeve OD. X3	Barrel length L	Total length "0" volume	Total length Full volume
PS-060037	20 ml	Ø29.0 [Ø1.14"]	Ø20.00 [Ø.7874"]	Ø33.35 [Ø1.312"]	133.35 [5.250"]	185.70 [7.312"]	250.80 [9.875"]
PS-060043	20 ml	Ø29.0 [Ø1.14"]	Ø20.00 [Ø.7874"]	Ø33.35 [Ø1.312"]	133.35 [5.250"]	185.70 [7.312"]	250.80 [9.875"]
PS-060039	50 ml	Ø33.5 [Ø1.32"]	Ø24.06 [Ø.9473"]	Ø37.85 [Ø1.490"]	179.40 [7.063"]	230.20 [9.063"]	338.10 [13.312"]
PS-060044	50 ml	Ø33.5 [Ø1.32"]	Ø24.06 [Ø.9473"]	Ø37.85 [Ø1.490"]	179.40 [7.063"]	230.20 [9.063"]	338.10 [13.312"]
PS-060040	100 ml	Ø39.4 [Ø1.55"]	Ø30.10 [Ø1.185"]	Ø43.80 [Ø1.725"]	206.40 [8.125"]	261.90 [10.312"]	403.20 [15.875"]
PS-060045	100 ml	Ø39.4 [Ø1.55"]	Ø30.10 [Ø1.185"]	Ø43.80 [Ø1.725"]	206.40 [8.125"]	261.90 [10.312"]	403.20 [15.875"]
PS-060041	250 ml	Ø55.9 [Ø2.20"]	Ø45.82 [Ø1.804"]	Ø60.35 [Ø2.375"]	231.80 [9.125"]	288.90 [11.375"]	441.30 [17.375"]
PS-060046	250 ml	Ø55.9 [Ø2.20"]	Ø45.82 [Ø1.804"]	Ø60.35 [Ø2.375"]	231.80 [9.125"]	288.90 [11.375"]	441.30 [17.375"]
PS-060042	500 ml	Ø66.0 [Ø2.60"]	Ø55.35 [Ø2.179"]	Ø71.10 [Ø2.800"]	298.45 [11.750"]	352.40 [13.875"]	560.40 [22.063"]
PS-060047	500 ml	Ø66.0 [Ø2.60"]	Ø55.35 [Ø2.179"]	Ø71.10 [Ø2.800"]	298.45 [11.750"]	352.40 [13.875"]	560.40 [22.063"]

Shown PS-060040

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.				Material MATERIAL		Dimensionen DIMENSIONS	
Aenderungen: MODIFICATIONS:				Projektion PROJECTION 		Allgemein Toleranzen GENERAL TOLERANCES SN EN 22768-1/ ISO 2768	
PART NO.:				Gezeichnet DRAWN 11.05.2010 F.P. Geprüft / Freigabe CHECKED/ RELEASE 11.05.2010 FCo		Masstab SCALE 0.8:1 CAD Nr. CAD NO. 506373.idw	
Magnum syringes General dimensions				Unterverz.: SUB DIREC.: KUNDE Zeichnungs Nr. / DRAWING NO.			
VICI AG International, CH-6214 Schenkon				3 - GE - 5063.73			