



SHOWN EU2



SHOWN EBU2L



Part#	Description	Bore		L		L2		L3		D		X1	X2	X3
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch			
EU2	External union 1/8" to 1/8"	Ø1.0	Ø.040"	-37.0	-1.456"							7/16" Hex		
EU2L	External union 1/8" to 1/8"	Ø2.0	Ø.080"	-37.0	-1.456"							7/16" Hex		
EU2T	External union 1/8" to 1/8"	Ø3.2	Ø.125"	-37.0	-1.456"							7/16" Hex		
EU4	External union 1/4" to 1/4"	Ø2.0	Ø.080"	-49.0	-1.930"							9/16" Hex		
EU4L	External union 1/4" to 1/4"	Ø4.6	Ø.180"	-49.0	-1.930"							9/16" Hex		
EU4T	External union 1/4" to 1/4"	Ø6.5	Ø.256"	-49.0	-1.930"							9/16" Hex		
EBU2L	Bulkhead external union 1/8" to 1/8"	Ø2.0	Ø.080"	-46.0	-1.811"	-12.2	~.480"	4.75	.187"	Ø8.0	Ø.3125"	7/16" Hex	7/16" Hex	1/2" Hex
EBU2T	Bulkhead external union 1/8" to 1/8"	Ø3.2	Ø.125"	-46.0	-1.811"	-12.2	~.480"	4.75	.187"	Ø8.0	Ø.3125"	7/16" Hex	7/16" Hex	1/2" Hex
EBU4	Bulkhead external union 1/4" to 1/4"	Ø2.0	Ø.080"	-54.6	-2.150"	-14.2	~.560"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	9/16" Hex	5/8" Hex
EBU4L	Bulkhead external union 1/4" to 1/4"	Ø4.6	Ø.180"	-54.6	-2.150"	-14.2	~.560"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	9/16" Hex	5/8" Hex
EBU4T	Bulkhead external union 1/4" to 1/4"	Ø6.5	Ø.256"	-54.6	-2.150"	-14.2	~.560"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	9/16" Hex	5/8" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / Modifications: a) 16.01.09 FP #568, added information

Material	Dimensions
Production	General Tolerances
EN 62278-1:ISO 2768	

Part No.:

**External unions including bulkhead versions**  
General dimensions

Scale	Order No.	07.03.2006	F.P.
5:1	Revision	16.01.2009	F.C.
	CAD No.	506072.idw	
Customer	KUNDE		
Drawing No.	Zeichnung Nr. / DRAWING NO.		

**VICI** VICI AG International, CH-6214 Schenkon **1 - GE - 5060.72 a**