


**SHOWN ZBU1**



Part#	Description	Bore		L		L2		L3		L4		D		X1	X2	X3
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch			
ZBU.5XC	1/32" union	Ø0.15	Ø.006"	~26.2	~1.030"	17.27	.680"	10.92	.430"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	3/16" Hex	1/2" Hex
ZBU.5	1/32" union	Ø0.25	Ø.010"	~26.5	~1.040"	17.53	.690"	11.18	.440"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	3/16" Hex	1/2" Hex
ZBU.5L	1/32" union	Ø0.50	Ø.020"	~27.0	~1.060"	18.03	.710"	11.68	.460"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	3/16" Hex	1/2" Hex
ZBU.5T	1/32" union	Ø0.89	Ø.035"	~25.7	~1.010"	16.76	.660"	10.41	.410"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	3/16" Hex	1/2" Hex
ZBU1XC	1/16" union	Ø0.15	Ø.006"	~36.8	~1.450"	21.21	.835"	14.86	.585"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU1C	1/16" union	Ø0.25	Ø.010"	~37.3	~1.470"	21.72	.855"	15.37	.605"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU1M	1/16" union	Ø0.50	Ø.020"	~37.8	~1.490"	22.23	.875"	15.88	.625"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU1	1/16" union	Ø0.75	Ø.030"	~38.3	~1.510"	27.73	.895"	16.38	.645"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU1L	1/16" union	Ø1.0	Ø.040"	~38.3	~1.510"	27.73	.895"	16.38	.645"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU1T	1/16" union	Ø1.65	Ø.065"	~36.4	~1.435"	20.83	.820"	14.48	.570"	3.18	.125"	Ø8	Ø.3125"	7/16" Hex	1/4" Hex	1/2" Hex
ZBU2	1/8" union	Ø0.75	Ø.030"	~50.1	~1.970"	31.12	1.225"	20.96	.825"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	3/8" Hex	5/8" Hex
ZBU2M	1/8" union	Ø1.0	Ø.040"	~50.1	~1.970"	31.12	1.225"	20.96	.825"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	3/8" Hex	5/8" Hex
ZBU2L	1/8" union	Ø2.0	Ø.080"	~50.1	~1.970"	31.12	1.225"	20.96	.825"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	3/8" Hex	5/8" Hex
ZBU2T	1/8" union	Ø3.25	Ø.128"	~48.2	~1.895"	29.21	1.150"	19.05	.750"	4.75	.187"	Ø11.1	Ø.4375"	9/16" Hex	3/8" Hex	5/8" Hex
ZBU4	1/4" union	Ø0.75	Ø.030"	~60.6	~2.385"	38.74	1.525"	24.77	.975"	6.35	.250"	Ø15.9	Ø.6250"	3/4" Hex	1/2" Hex	7/8" Hex
ZBU4L	1/4" union	Ø4.6	Ø.180"	~60.6	~2.385"	38.74	1.525"	24.77	.975"	6.35	.250"	Ø15.9	Ø.6250"	3/4" Hex	1/2" Hex	7/8" Hex
ZBU4T	1/4" union	Ø6.48	Ø.255"	~58.7	~2.310"	36.83	1.450"	22.86	.900"	6.35	.250"	Ø15.9	Ø.6250"	3/4" Hex	1/2" Hex	7/8" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / Modifications: a) 15.01.09 FP #568, added information

Material	Dimensions
Production	General Tolerances

Part No. **Bulkhead internal unions**  
**Stainless steel**  
 General dimensions

Scale: 5:1  
 Date: 27.03.2006  
 CAD No.: 15.01.2009  
 Drawing: 506079.idw  
 Customer: KUNDE

VICI® VICI AG International, CH-6214 Schenkon  
 1 - GE - 5060.79 a